

Moxi™ at-a-glance

Styled for:



Performance profile

	premium	e16	e8	essential
Range of loss	Mild to severe	Mild to severe	Mild to severe	Mild to severe
Channels	16	16	8	4
Bands	16	16	8	8
Processing types	WDRC	WDRC Linear Limiting	WDRC Linear Limiting	WDRC Linear Limiting

Features

	AutoPro4™ (4 sound destinations)	AutoPro2™ (2 sound destinations)	AutoMic Automatic directional system (fixed or adaptive)	
Automatic programs				
Manual listening programs	Up to 3	Up to 3	Up to 3	Up to 3
OnBoard™ control (VC or program control)	✓	✓	✓	✓
Directional microphone system	Multiband adaptive directionality	Adaptive directionality	Adaptive directionality	Fixed directional
AntiShock™	✓	✓	✓	✓
Speech enhancement	Speech enhancement LD	Speech enhancement LD	✓	
Phase canceller	✓	✓	✓	✓
Data logging	including VC changes	including VC changes	✓	✓
MyMusic™	✓			
Noise reduction	✓	✓	✓	✓
Wind noise manager	✓	✓	✓	✓
Quiet mode expansion	✓	✓	✓	✓
Easy-t	✓	✓	✓	✓
Telecoil	✓	✓	✓	✓
Program beep indicator	✓	✓	✓	✓
Low battery warning	✓	✓	✓	✓
Multiple choices of instant fit	✓	✓	✓	✓



Four products. Four price points.
One CRT product family.